ČZ - INFORMACE PRO UŽIVATELE

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EN - USER INFORMATION

FLOWER - powder free nitrile examination gloves blue, black, non-sterile. Available size: XS, S, M, L, XL.

1) Medical Device Regulation (MDR) (EU) 2017/745

- This product is classified under Class I Medical Device per Rule 1 and Rule 5 of Annex IX, meets the provisions of the Medical Device Regulation (MDR) (EU) 2017/745.
- This product conforms with European Standards EN 455-1:2020/A1:2022, EN 455-2:2015, EN 455-3:2015.

2) PPE EU Type-Examination

- This product is classified as Category III Personal Protective Equipment (PPE) according to PPE Regulation (EU) 2016/425 and has been shown to comply with this Regulation through the Harmonised European Standards EN ISO 21420, EN ISO 374-1:2016+A1:2018/Type B, EN ISO 374-4:2019 a EN ISO 374-5:2016
- Notified Body responsible for certification is SATRA Technology Europe Limited (2777), Bracetown Business Park, Clonee, Dublin 15, D15 YN2P, Ireland.
- Notify Body responsible for internal production control plus supervised product checks at random intervals (Module C2) is SATRA Technology Europe Limited (2777), Bracetown Business Park, Clonee, Dublin 15, D15 YN2P, Ireland.
- The EU Declaration of Conformity is accessible at www.batist.com

- Micro Organism Hazards Pictogram: EN ISO 374-5:2016 Protect against Bacteria, Fungi and Virus. No pathogens have been detected through the specimen and the following pictogram is applied.

- Chemical Hazards Pictogram: EN 16523-1:2015; Additional information on chemical resistance obtainable from manufacturer.

EN ISO 374-1:2016+A1:2018/Type B permeation levels are based on breakthrough times as follows:

Permeation Performance Level	1	2	3	4	5	6
Measured breakthrough time (min)	>10	>30	>60	>120	>240	>480

This product complies with Type B requirements and the following pictogram shall be used with reference to clause 6.2 of ISO 374-1.

- Performance and Limitation of Use
- This product has been tested in accordance with EN ISO 374-5:2016. Protection against bacteria and fungi - Pass Protection against viruses - Pass

- Gloves had been tested in accordance with EN16523-1:2015+A1:2018 a EN 16523-2:A1:2018 resistance to permeation by chemicals and achieved the following performance levels:

Chemicals	Performance Level
40% Sodium Hydroxide (K)	6
37% Formaldehyde (T)	6
30% Hydrogen peroxide (P)	6

- This information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals.
- The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only (except in cases where the glove is equal to or over 400 mm² - where the cuff is tested also) and relates only to the chemical tested. It can be different if the chemical used is in a mixture.
- It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation.
- When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves.
- Before usage, inspect the gloves for any defect or imperfections.
- This product had been tested in accordance with EN 374-4:2019 and achieved the following degradation results:

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